



X-BOLT®

The Hip Fracture Fixation Company



X-BOLT Hip System (XHS™)
Surgical Technique





X-BOLT Mini-Plate
Surgical Technique

X-BOLT® SIMPLE AND SECURE


The X-BOLT® provides strong femoral head fixation and excellent rotational stability. The X-BOLT® comprehensive systems replicate traditional techniques. The simple, concise and user-friendly instrumentation ensures fast and dependable hip fracture surgery.

The user-friendly features allow surgeons the freedom to concentrate on the important steps of:

- ✓ **Fracture Reduction**
- ✓ **Guidewire Position**



X-BOLT® HOW IT WORKS




The X-BOLT® is an expanding bolt akin to a Chinese lantern with a central drive shaft. The opposing threads compress the expandable section from both ends to expand the wings perpendicularly to the shaft, without spinning, pushing or pulling the femoral head. Maximal expansion of the four wings is always easily seen on AP and lateral fluoroscopy.

- ✓ **Expanded (or retracted) with standard screwdriver**
- ✓ **Compacts surrounding cancellous bone**
- ✓ **Maximises resources within femoral head**
- ✓ **Single centre-centre guidewire shot**
- ✓ **No need for de-rotation screw**
- ✓ **Easily reversible and removable**

XHS™ DYNAMIC HIP PLATING

1. **SET-UP.** Follow all steps on fluoroscopy and **check rotation or lateral view if any step is troublesome.** Reduce fracture on fracture table. Anatomical reduction is key to successful outcome. Standard approach through vastus lateralis to lateral femur.
2. **GUIDEWIRE.** Using angle guide, place femoral head guidewire centre-centre to within 2-3 mm of joint line (optimum tip-apex point). Remove angle guide. **In unstable fracture patterns, consider placing a second guidewire** to skewer and anchor reduced fracture position, as primary guidewire will lose its grip when over-reamed.
3. **MEASURE.** Measure depth of laser mark on guidewire shaft with guidewire measuring ruler. **Round down to nearest 5mm** for reamer and X-Bolt length. Note X-Bolt shortens by 2mm from base end on full expansion. Make allowance for any post-operative fracture collapse.
4. **REAM.** Set reamer to measured length. Insert over guidewire. Ream to tip-apex point, also creating barrel channel and lateral counter-sink. Follow closely under fluoroscopy, ensuring guidewire does not advance inadvertently. **Ream fully deep to tip-apex point** (reset longer depth if necessary), as final drilling depth dictates X-Bolt final position. Remove reamer and primary guidewire.
5. **BONE CRUSHER.** Insert bone crusher fully to tip-apex point. Starting with handle facing at 6 o'clock position, deploy at 'hourly' intervals, compacting the cancellous bone, through an arc of 180°, to the 12 o'clock position. This creates a diamond-shaped starter cavity that the X-Bolt will expand into. **Note: It is always necessary to complete this step regardless of the nature/density of bone.**
6. **INSERT X-BOLT AND PLATE.** On side table, slide X-Bolt into plate barrel. Insert screwdriver into base of X-Bolt for use as an insertion instrument; the screw holding sleeve prevents plate from sliding off. Insert X-Bolt and plate into reamed channel. Remove the screw holding sleeve and advance X-Bolt fully to tip-apex point. If necessary, seat plate onto femur using plastic impactor. Remove the second guidewire if it has been used.
7. **STABILISE PLATE.** Drill, measure and place most proximal cortical screw through the plate, ensuring correct plate rotation on femur. The stabilised plate counters X-Bolt rotation when deploying.
8. **EXPAND X-BOLT.** Follow expansion on fluoroscopy. Expand by clockwise screwdriver rotation until a stop is felt, or to 4.5Nm torque screwdriver limit, or by surgeon choice. If initial resistance to expansion is felt, do not force; remove plate and repeat bone crusher step. **Note: X-Bolt must be discarded and replaced if expansion has been reversed, or if the X-Bolt is damaged in any way.**
9. **(Optional) FRACTURE COMPRESSION.** Remove leg traction. Insert compression rod (note: anticlockwise threads) into the base of X-Bolt. Manually retract rod and thus X-Bolt, using plastic impactor as counter force, to achieve desired compression.
10. **FINAL STEPS.** Insert remaining cortical screws and check final position. Wound closure in layers, post-operative management and weight-bearing as per surgeon's instructions.

XDG-135
*XNI-010

XMG-001

XRT-090

XBC-001

XHS-035
XSH-001

XPI-300

XCG-035
*XDS-032
XNI-030

XHS-035

XSR-300
XPI-300

XHS-035

X-BOLT® MINI-PLATE

Part Code



- 1. INDICATIONS.** The X-Bolt Mini-Plate should be reserved for patients with stable intracapsular fractures (Garden I & II).
- 2. SET-UP.** Follow all steps on AP & lateral fluoroscopy. Use a mini-open approach through or below vastus lateralis to lateral femur.
- 3. GUIDEWIRE.** Freehand starting point: opposite superior aspect of lesser trochanter. Place femoral head guidewire centre-centre in femoral head to within 2-3mm of joint line (optimum tip-apex point). The green & gold drill guide can be used as a soft tissue protector.
- 4. MEASURE.** Measure depth of laser mark on guidewire shaft with guidewire measuring ruler. Note: X-Bolt shortens by 2mm from base end on full expansion and to ensure that X-Bolt does not disengage in situ with mini-plate, **add 5-10mm for X-Bolt length.**
- 5. REAM.** Remove barrel part of triple reamer. Ream over guidewire to tip-apex point, to within 2-3mm of joint line (optimum tip-apex point). **Follow closely under fluoroscopy** to ensure guidewire or reamer does not advance inadvertently. Remove reamer and guidewire.
- 6. BONE CRUSHER.** Insert bone crusher fully to tip-apex point. Starting with handle facing at 6 o'clock position, deploy at 'hourly' intervals, compacting the cancellous bone, through an arc of 180°, to the 12 o'clock position. This creates a diamond-shaped starter cavity that the X-Bolt will expand into. **Note: It is always necessary to complete this step regardless of the nature/density of bone.**
- 7. INSERT X-BOLT AND MINI-PLATE.** On side table, slide X-Bolt into mini-plate. Insert screwdriver into base of X-Bolt for use as an insertion instrument; the screw holding sleeve prevents mini-plate from sliding off. Insert X-Bolt into reamed channel. Remove the screw holding sleeve and advance X-Bolt fully to tip-apex point.
- 8. LOCKING SCREW.** Insert threaded drill guide into hole in plate. Drill via guide using Ø4.0mm drill bit. This should not penetrate medial cortex. Remove drill guide. Using screwdriver, place unicortical 30mm x Ø5.0mm locking screw into plate. This stabilises the mini-plate to counter rotation of X-Bolt when deploying.
- 9. EXPAND X-BOLT.** Follow expansion on fluoroscopy. Expand by clockwise screwdriver rotation until a stop is felt, or to 4.5Nm torque screwdriver limit, or by surgeon choice. If initial resistance to expansion is felt, do not force; remove mini-plate and X-Bolt and repeat bone crusher step. **Note: X-Bolt must be discarded and replaced if expansion has been reversed, or if the X-Bolt is damaged in any way.**
- 10. CLOSURE.** Wound closure, post-operative management and weight-bearing as per surgeon's instructions.

*XNI-010
XCG-035

XMG-001

XRT-090

XBC-001

*XMP-125
XHS-035
XSH-001

*XNI-029

XHS-035

TIPS FOR SUCCESS

1. **ANATOMICAL REDUCTION.** Avoid varus; check medial cortex congruity. Lateral view on fluoroscopy vital to identify and correct sag and to check rotational alignment.
2. **REAMING TO TIP-APEX POINT.** Vitaly important to maintain a deep tip-apex point to further reduce the risk of cut-out.
 - (a) Reamer depth dictates X-Bolt final position
 - (b) Bone crusher must be inserted fully to tip-apex point
3. **BONE CRUSHER.** Deploy at multiple intervals to compact cancellous bone, creating a 18mm span diamond-shaped starter cavity. No bone is removed, it is compacted radially. It further improves rotational stability by creating a peripheral zone of compacted bone. The X-Bolt expands to a maximum 24mm span.



Part Code

REMOVAL INSTRUCTIONS

- | | |
|---------------------|--|
| XHS-035 | 1. SCREWDRIVER. Reverse X-Bolt expansion with anticlockwise rotation of screwdriver inserted into base of X-Bolt. The X-Bolt is able to crush new bone formed within expanded wings. Remove cortical screws and plate when X-Bolt wings fully retracted. |
| XSR-300 | 2. COUPLING SCREW. Insert coupling screw rod into X-Bolt, noting anticlockwise threads. Remove X-Bolt. |
| _____ | |
| XHS-035 | 1. BROKEN REMOVAL. In the rare event of X-Bolt breaking at the wing elbows, continue anticlockwise turns to disengage drive screw from deep parachute fragment. |
| XSR-300 | 2. REMOVE SHAFT OF X-BOLT. Insert holder coupling screw rod into base of X-Bolt, noting the anticlockwise threads. Remove X-Bolt shaft and drive screw, along with cortical screws and plate, leaving the parachute fragment in-situ. |
| | 3. REMOVE PARACHUTE FRAGMENT WITH FEMORAL HEAD. Parachute fragment may be left in-situ or removed with femoral head if converting to hemiarthroplasty or total hip replacement. |
| *XNI-038 XNI-037 | 4. PARACHUTE GRASPER. If necessary to remove parachute fragment, insert threaded grasper guidewire into parachute fragment. Slide parachute grasper over guidewire, under fluoroscopy, and flatten parachute wings. The parachute fragment and guidewire are removed together by manually pulling backwards on the guidewire. |

STANDARD REMOVAL



BROKEN REMOVAL



(* single
packed
sterile)

ORDERING INFORMATION

X-BOLT®
stainless steel
Ø 9.0mm

| Code | NHS NPC | Description |
|---------|---------|--------------|
| XBB-080 | FJY7494 | 80mm X-Bolt |
| XBB-085 | FJY7495 | 85mm X-Bolt |
| XBB-090 | FJY7496 | 90mm X-Bolt |
| XBB-095 | FJY7497 | 95mm X-Bolt |
| XBB-100 | FJY7498 | 100mm X-Bolt |
| XBB-105 | FJY7499 | 105mm X-Bolt |
| XBB-110 | FJY7500 | 110mm X-Bolt |
| XBB-115 | FJY7501 | 115mm X-Bolt |
| XBB-120 | FJY7502 | 120mm X-Bolt |
| XBB-125 | FJY7503 | 125mm X-Bolt |

XHS™ Barrel-Plate
stainless steel
135-degree

| | | |
|----------|---------|----------------------------------|
| XBP-1352 | FJY7504 | 2-hole 135° Dynamic Barrel-Plate |
| XBP-1354 | FJY7505 | 4-hole 135° Dynamic Barrel-Plate |
| XBP-1355 | FJY7506 | 5-hole 135° Dynamic Barrel-Plate |
| XBP-1356 | FJY7507 | 6-hole 135° Dynamic Barrel-Plate |

X-BOLT® Mini-Plate
stainless steel
125-degree

| | | |
|---------|---------|--|
| XMP-125 | FJY7509 | Mini-Plate, 30mm x Ø5.0mm locking screw, locking screw drill guide |
|---------|---------|--|

Cortical Screws
stainless steel
Ø 4.5mm, self-tapping

| | | |
|---------|---------|--|
| XSS-030 | FJY7535 | 30mm x Ø4.5mm cortical screw, self-tapping |
| XSS-032 | FJY7536 | 32mm x Ø4.5mm cortical screw, self-tapping |
| XSS-034 | FJY7537 | 34mm x Ø4.5mm cortical screw, self-tapping |
| XSS-036 | FJY7538 | 36mm x Ø4.5mm cortical screw, self-tapping |
| XSS-038 | FJY7539 | 38mm x Ø4.5mm cortical screw, self-tapping |
| XSS-040 | FJY7540 | 40mm x Ø4.5mm cortical screw, self-tapping |
| XSS-042 | FJY7541 | 42mm x Ø4.5mm cortical screw, self-tapping |
| XSS-044 | FJY7542 | 44mm x Ø4.5mm cortical screw, self-tapping |
| XSS-046 | FJY7543 | 46mm x Ø4.5mm cortical screw, self-tapping |
| XSS-048 | FJY7544 | 48mm x Ø4.5mm cortical screw, self-tapping |

Femoral Head Guidewire

| | | |
|---------|---------|------------------------------------|
| XNI-010 | FJY7519 | Ø 3.2mm x 385mm threaded guidewire |
|---------|---------|------------------------------------|


Drill bit for XHS™ cortical screws

| | | |
|---------|---------|--|
| XDS-032 | FJY7508 | Ø 3.2mm x 145mm drill bit, quick connect |
|---------|---------|--|


Drill bit for mini-plate locking screw

| | | |
|---------|---------|--|
| XNI-029 | FJY7521 | Ø 4.0mm x 150mm pilot drill bit, quick connect |
|---------|---------|--|

INSTRUMENTATION (SET P)



| Part Code | Description |
|------------------|--|
| XRT-090 | Triple Reamer |
| XCG-035 | Drill Guide, Green & Gold, Ø 3.2mm |
| XDG-135 | 135-degree drill guide |
| XHP-035 | Power driver, 3.5mm-hex, quick coupling |
| XSH-001 | Screw Holding Sleeve |
| XMG-001 | Guidewire Measuring Ruler |
| XBC-001 | Bone Crusher |
| XPI-300 | Plastic Impactor |
| XNI-030 | Depth Gauge for cortical screws |
| XHS-035 | Screwdriver, 3.5mm-hex, 4.5Nm torque limit |
| XSR-300 | Compression Screw Rod |



Set Z (Removal Set)

| Part Code | Description |
|------------------|---------------------------------|
| XNI-033 | X-Bolt extractor shaft |
| XNI-034 | Nail extractor shaft |
| XNI-035 | Slap hammer |
| XNI-036 | Crochet hook removal tool |
| XNI-037 | Parachute grasper |
| XNI-038 | Guidewire for parachute grasper |



Your local contact:

X-BOLT DIRECT LTD
+44 1172 300632

+44 7585 639414 (mobile)
+44 7825 372750 (mobile)
+44 7469 879150 (mobile)

For further information, or to order, please contact:



X-Bolt Orthopaedics
Unit 5, Northwood Court
Santry, Dublin 9, Ireland



+353 1 8456011



info@x-bolt.com



www.x-bolt.com

The information presented in this brochure is intended as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific X-Bolt products. Always refer to the package insert, product label and instructions for use before using any X-Bolt product. Surgeons must always rely on their own clinical judgement, training and expertise when deciding which products and techniques to use with their patients.

X-Bolt® is a registered trademark of X-Bolt Orthopaedics

SOTA Orthopaedics Limited, trading as X-Bolt Orthopaedics, is registered in Ireland, number 439651

European Patent No. EP 2175790 and US Patent Application No. 12/667,513

SURGT-XBOLT-003, Revision 1

© 2016 X-Bolt Orthopaedics, all rights reserved.

