



IM Hip Nailing Surgical Technique



SURGICAL TECHNIQUE

- XNI-028
- 1. NAIL SIZE, LEFT OR RIGHT** All nails are sided left or right. For long nails, measure nail length by the distance from the greater trochanter to the superior aspect of patella with long nail ruler, under fluoroscopy.
- XNI-019
XNI-013
XNI-007
- 2. JIG ASSEMBLY.** Mount nail to jig using M-10 connection bolt and flexi-drive. Remove flexi-drive. Manually check alignment is correct with drill sleeves.
- *XNI-010
XNI-004
XNI-005
XNI-002
*XNI-001
- 3. REDUCTION.** Anatomically reduce fracture prior to insertion of nail. Follow all steps on fluoroscopy.
- 4. ENTRY POINT. Medial to the tip of greater trochanter.** Make a small proximal skin incision in line with expected path. Use Ø3.2mm femoral head guidewire and proximal sleeves, or cannulated awl. Pass long guidewire down intramedullary canal into distal fragment.
- XNI-006
XNI-004
XNI-005
- 5. REAM.** Use conical reamer to ream proximal femur, with soft tissues protected with proximal sleeves. If distal reaming needed, use flexible reamers in Ø0.5mm increments to Ø12.0mm (Ø2.0mm greater than nail diameter).
- 6. NAIL INSERTION.** Insert nail coupled to jig over guidewire (having removed flexi-drive) and advance by hand. **DO NOT HIT alignment jig.** Remove guidewire when nail properly in distal fragment. Gentle final impaction may be performed only if absolutely necessary, onto and via flexi-drive re-inserted into jig, after long guidewire is removed.
- XNI-008
XNI-009
*XNI-010
- 7. FEMORAL HEAD GUIDEWIRE.** Through the jig and via stab incision, insert inner and outer X-Bolt sleeves. Place Ø3.2mm femoral head guidewire centre-centre in femoral head within 2-3mm of joint line cortex; optimum tip-apex point, checking position on AP and lateral.
- XNI-011
- 8. MEASURE.** Use guidewire ruler to measure guidewire protrusion from inner sleeve. Note: X-Bolt shortens by 2-3mm from base-end on expansion. Do not bury X-Bolt deep to lateral cortex that could prevent dynamic sliding. **Add approximately +5mm for X-Bolt length.**
- XNI-012
- 9. REAM.** Remove inner sleeve. Set Ø9.1mm reamer depth, at base end of stopper. Ream over guidewire to create X-Bolt channel to within 2-3mm of joint line (tip-apex point). Leave reamer temporarily in situ, (but detach from power driver) to prevent alignment jig from sagging/rotation from gravity.
- XNI-015
XNI-016
*XNI-017
XNI-018
- 10. DISTAL LOCKING. Short nail:** Insert green and gold sleeves through appropriate (left/right) aperture, choosing dynamic or static distal locking option. Use Ø4.0mm x 305mm drill bit via inner (gold) sleeve. Remove inner (gold) sleeve. Confirm screw length using measuring hook. Insert appropriate Ø5.0mm distal locking screw. **Long Nail:** Use Ø4.0mm x 150mm drill bit for freehand distal locking, in dynamic or static configurations. Screw length is measured with standard depth gauge.
- *XNI-029
XNI-030





11. **INSERT BONE CRUSHER.** Remove reamer and guidewire. This leaves a clear bone channel in femoral head. Insert bone crusher fully deep, to within 2-3mm of joint line (tip-apex point).

XNI-014



12. **DEPLOY BONE CRUSHER.** Deploy liberally at multiple intervals to compact cancellous bone through an arc of 180°, thus creating a small diamond-shaped 360° starter cavity (span 18mm) in femoral head. The X-Bolt subsequently can expand to 24mm span. **Note: It is always necessary to complete this step, regardless of the nature/density of bone.**

13. **INSERT X-BOLT.** Mount the sized X-Bolt onto T-handle screwdriver and advance X-Bolt into position in femoral head with flat aspect facing superiorly.

XNI-022



14. **INSERT FLEXI-DRIVE AND FLEXIBLE SCREWDRIVER.** Insert flexi-drive into hollow metal proximal end of jig. Insert flexible screwdriver through proximal end of flexi-drive to engage pre-loaded set screw within the nail.

XNI-007

XNI-023



15. **ENGAGE SET-SCREW.** Gently advance set screw (light finger grip on screwdriver) to just lightly kiss onto the superior flat aspect of X-Bolt. The set screw prevents rotation, but should not prevent dynamic sliding, which can be confirmed by gently rocking X-Bolt via T-handle screwdriver prior to expansion. **Option:** The set screw may be fully tightened down as a final step following X-Bolt expansion if a locked construct is desired.

16. **EXPAND X-BOLT.** Insert T-handle torque-limiting screwdriver into base of X-Bolt. Standard clockwise turns expand the X-Bolt until a stop is felt, or the 4.5Nm torque limit is reached. Do not use excessive force. If inadvertent resistance is felt, remove the X-Bolt and repeat bone crusher step. Reversal, if necessary, is performed by anticlockwise screwdriver rotation. **NB: X-Bolt must be discarded and replaced if full expansion has been reversed, or if X-Bolt is damaged in any way.**



17. **END-CAP.** Remove jig by de-coupling connector bolt via flexible coupler, using screwdriver shaft through horizontal apertures if greater leverage is required. Reinsert Ø2.4mm x 900mm long guidewire into top of nail. Place appropriate +2mm or +5mm end-cap over long guidewire into top of nail using cannulated 4.0mm-hex screwdriver. Remove long guidewire.

*XNI-001

XNI-024



18. **CLOSURE.** Wound closure and post-operative management as per surgeon's instructions.


(* = single
packed
sterile)

X-BOLT® : ACHIEVING GREAT RESULTS

The X-Bolt provides strong femoral head fixation and excellent rotational stability, whilst the nail (14.5mm proximal diameter) causes less disruption to the greater trochanter and hip abductors. The novel jig and flexi-drive, pre-loaded set-screw and cannulated end-cap make for simple and fast surgery, without fiddle-factor.

These user-friendly features allow surgeons the ability and freedom to concentrate on the important surgical aspects of:

- ✓ **Fracture Reduction**
- ✓ **Nail Entry Point**
- ✓ **Femoral Head Guidewire Position**



TROUBLESHOOTING

Problem	Possible Cause	Solution
Jig not connecting properly to nail	Set-screw too high in nail	Advance set-screw with T-Handle or flexible screwdriver (3.5mm hex)
Difficulty passing long guidewire or nail across fracture site	Inadequate fracture reduction	Accurate reduction should be performed, with a low threshold for open reduction or minimally invasive fracture reduction techniques
Femoral head guidewire will not pass through nail	Long guidewire still in situ	Remove long guidewire
Ø9.1mm reamer will not pass through nail	Set-screw blocking progress	Insert flexi-drive and flexible screwdriver into proximal end of jig and retract the set-screw superiorly
Ø9.1mm reamer will not pass through nail	Jig sagged due to gravity and lever-arm ('Jig sag')	Leaving the Ø9.1mm reamer temporarily in-situ will help stabilise jig, as does the distal locking screw when implanted
Ø4.0mm drill-bit won't pass through nail	Not aligned with the hole	Tighten jig to nail Ensure not leaning on jig and/or bending jig
Cannot implant distal locking screw after drilling	Jig sag	Leave Ø9.1mm reamer in place until distal locking is complete
Bone crusher unable to advance to deep tip-apex point	Too shallow reaming with Ø9.1mm reamer	Re-insert guidewire and ream again to deep tip-apex point with Ø9.1mm reamer
Cannot implant X-Bolt after reaming	Jig sag	Distal lock before implanting X-Bolt
X-Bolt spins when attempting to deploy	Set-screw not engaged on flat aspect of X-Bolt shaft	Rotate X-Bolt until flat aspect faces superiorly and then gently advance set screw (light finger grip on screwdriver) to just lightly kiss onto the superior flat aspect of X-Bolt

REMOVAL INSTRUCTIONS

Part Code

STANDARD REMOVAL



1. **STANDARD REMOVAL.** Retract X-Bolt expansion with anticlockwise rotation of screwdriver, inserted into base of X-Bolt. The retraction of the X-Bolt will crush new bone formed within expanded wings.
2. **REMOVE END-CAP.** When wings fully retracted, insert guidewire into cannulated end-cap. Remove end-cap using 4.0mm-hex cannulated screwdriver.
3. **LOOSEN SET-SCREW.** Insert 3.5mm-hex screwdriver and loosen set screw sufficiently to enable free passage of X-Bolt through nail aperture.
4. **REMOVE X-BOLT.** Insert removal rod into X-Bolt, noting the anticlockwise (left hand) threads. Remove X-Bolt using removal rod, and slap hammer extension, if necessary.
5. **REMOVE NAIL.** Remove distal locking screw with standard 3.5mm hex screwdriver. Remove nail using removal rod and slap hammer, if necessary.

XNI-022

XNI-024

XNI-033
XNI-035

XNI-034
XNI-035

BROKEN REMOVAL



1. **BROKEN REMOVAL.** In the rare event of X-Bolt breaking at the wing elbows, continue anticlockwise turns to disengage drive screw from deep parachute fragment.
2. **END-CAP AND SET-SCREW.** Remove end-cap and loosen set-screw, as described above.
3. **REMOVE SHAFT OF X-BOLT.** Insert removal rod into base of X-Bolt, noting the anticlockwise (left hand) threads. Remove X-Bolt shaft and drive screw, leaving the parachute fragment in-situ.
4. **REMOVE NAIL.** Remove distal locking screw and nail, as described above.
5. **PARACHUTE FRAGMENT.** If removal of parachute fragment is desired (it may be left in-situ or removed with femoral head if converting to arthroplasty), place grasper guidewire into right-hand threads of the parachute fragment until a stop is felt.
6. **PARACHUTE GRASPER.** Slide grasper over guidewire, under fluoroscopy and flatten parachute wings in orthogonal planes. The parachute fragment is removed by manually pulling backwards on the guidewire.

XNI-022

XNI-024

XNI-033
XNI-035

XNI-034
XNI-035

*XNI-038

XNI-037

ORDERING INFORMATION

X-Bolt®
stainless steel
Ø 9.0mm

Code	Description
XBB-080	80mm X-Bolt
XBB-085	85mm X-Bolt
XBB-090	90mm X-Bolt
XBB-095	95mm X-Bolt
XBB-100	100mm X-Bolt
XBB-105	105mm X-Bolt
XBB-110	110mm X-Bolt
XBB-115	115mm X-Bolt
XBB-120	120mm X-Bolt
XBB-125	125mm X-Bolt

IM Hip Nail
stainless steel
130-degree
4° valgus
14.5mm proximal diameter
10.0mm distal diameter
2.0m radius bend for long nails
All sided Left or Right
Pre-loaded set-screw

XNN-191	195mm -Short- IM Hip Nail (left)
XNN-192	195mm -Short- IM Hip Nail (right)
XNN-341	340mm X-Bolt Nail (left)
XNN-342	340mm X-Bolt Nail (right)
XNN-361	360mm X-Bolt Nail (left)
XNN-362	360mm X-Bolt Nail (right)
XNN-381	380mm X-Bolt Nail (left)
XNN-382	380mm X-Bolt Nail (right)
XNN-401	400mm X-Bolt Nail (left)
XNN-402	400mm X-Bolt Nail (right)
XNN-421	420mm X-Bolt Nail (left)
XNN-422	420mm X-Bolt Nail (right)

Nail End-Cap
stainless steel
Cannulated 4.0mm-hex socket


XNE-002	End-cap +2.0mm
XNE-005	End-cap +5.0mm

Distal Locking Screw
stainless steel
Ø 5.0mm, self-tapping

XND-030	30mm x Ø5.0mm Distal locking screw
XND-035	35mm x Ø5.0mm Distal locking screw
XND-040	40mm x Ø5.0mm Distal locking screw
XND-045	45mm x Ø5.0mm Distal locking screw
XND-050	50mm x Ø5.0mm Distal locking screw
XND-055	55mm x Ø5.0mm Distal locking screw

Single Sterile Packaged Consumables

XNI-001	Ø 2.4mm x 900mm long guidewire
XNI-010	Ø 3.2mm x 385mm femoral head guidewire
XNI-017	Ø 4.0mm x 305mm pilot drill bit, quick connect
XNI-029	Ø 4.0mm x 150mm pilot drill bit, quick connect
XNI-038	Parachute grasper guidewire




Set X

Description


Part Code

- Curved Awl Cannulated XNI-002
- T-handle quick-fit Jacobs chuck XNI-003
- Outer proximal reamer sleeve XNI-004
- Inner proximal reamer sleeve XNI-005
- Rigid conical reamer XNI-006
- Flexi-drive for M-10 bolt XNI-007
- M-10 connector bolt XNI-013
- Alignment jig XNI-019
- Long nail measuring ruler XNI-028



Set Y

- Outer drill sleeve for X-Bolt XNI-008
- Inner drill sleeve for X-Bolt XNI-009
- Guidewire ruler for X-Bolt XNI-011
- Ø9.1mm X-Bolt drill-reamer XNI-012
- Bone crusher XNI-014
- Outer distal locking sleeve XNI-015
- Inner distal locking sleeve XNI-016
- Depth hook for locking screw XNI-018
- Compression screw rod XNI-021
- 3.5mm-hex screwdriver (T-handle) (4.5Nm torque limit) XNI-022
- 3.5mm-hex screwdriver (flexible) XNI-023
- 4.0mm-hex screwdriver (cannulated) XNI-024
- Depth gauge for cortical screw XNI-030



Set Z (Removal set)

- X-Bolt extractor shaft XNI-033
- Nail extractor shaft XNI-034
- Slap hammer XNI-035
- Crochet hook removal tool XNI-036
- Parachute grasper XNI-037
- Guidewire for parachute grasper XNI-038




Your local contact:

**X-BOLT DIRECT LTD
BRISTOL & BATH SCIENCE PARK
EMERSON'S GREEN, BRISTOL
BS16 7FR**

**+44 1172 300632 (office)
+44 78 25372750 (mobile)**

For further information, or to order, please contact:

 X-Bolt Orthopaedics
Unit 5, Northwood Court
Santry, Dublin 9, Ireland

 **+353 1 8456011**
 **info@x-bolt.com**
 **www.x-bolt.com**

The information presented in this brochure is intended as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific X-Bolt products. Always refer to the package insert, product label and instructions for use before using any X-Bolt product. Surgeons must always rely on their own clinical judgement, training and expertise when deciding which products and techniques to use with their patients.

X-Bolt® is a registered trademark of X-Bolt Orthopaedics

SOTA Orthopaedics Limited, trading as X-Bolt Orthopaedics, is registered in Ireland, number 439651

European Patent No. EP 2175790 and US Patent Application Nos. 12/667,513 & 13/177,411

SURGT-XBOLT-002, Revision 2

© 2014 X-Bolt Orthopaedics, all rights reserved.

